


## alpha-Adenosine

## Chemical Properties

CAS No. :	5682-25-7
Formula:	C <sub>10</sub> H <sub>13</sub> N <sub>5</sub> O <sub>4</sub>
Molecular Weight:	267.24
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	alpha-Adenosine is a useful organic compound for research related to life sciences. The catalog number is TNU1637 and the CAS number is 5682-25-7.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.742 mL	18.7098 mL	37.4195 mL
5 mM	0.7484 mL	3.742 mL	7.4839 mL
10 mM	0.3742 mL	1.871 mL	3.742 mL
50 mM	0.0748 mL	0.3742 mL	0.7484 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481