

Compound N016-0008

Chemical Properties


CAS No. :

Formula: C₅₄H₈₆O₁₆

Molecular Weight: 991.266

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N016-0008 is a useful organic compound for research related to life sciences and the catalog number is T131352.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0088 mL	5.044 mL	10.0881 mL
5 mM	0.2018 mL	1.0088 mL	2.0176 mL
10 mM	0.1009 mL	0.5044 mL	1.0088 mL
50 mM	0.0202 mL	0.1009 mL	0.2018 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481