


## Thianaphthene-2-carboxylic acid

## Chemical Properties

CAS No. : 6314-28-9  
Formula: C<sub>9</sub>H<sub>6</sub>O<sub>2</sub>S  
Molecular Weight: 178.21  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Thianaphthene-2-carboxylic acid, with CAS No. 6314-28-9, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Thianaphthene-2-carboxylic acid provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.6114 mL	28.0568 mL	56.1136 mL
5 mM	1.1223 mL	5.6114 mL	11.2227 mL
10 mM	0.5611 mL	2.8057 mL	5.6114 mL
50 mM	0.1122 mL	0.5611 mL	1.1223 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481