

## Compound T127923

## Chemical Properties

CAS No. :

Formula: C<sub>20</sub>H<sub>28</sub>O<sub>7</sub>

Molecular Weight: 380.44

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	(2E,4E)-6-((3-hydroxy-3,5-dimethyl-6-[(4E)-4-methylhex-4-en-2-yl]-2-oxooxan-4-yl)oxy)-6-oxohexa-2,4-dienoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6285 mL	13.1427 mL	26.2854 mL
5 mM	0.5257 mL	2.6285 mL	5.2571 mL
10 mM	0.2629 mL	1.3143 mL	2.6285 mL
50 mM	0.0526 mL	0.2629 mL	0.5257 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481