


2,6-Dichloroisonicotinic acid

Chemical Properties

CAS No. :	5398-44-7
Formula:	C ₆ H ₃ Cl ₂ NO ₂
Molecular Weight:	192
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	2,6-Dichloroisonicotinic acid, with CAS No. 5398-44-7, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 2,6-Dichloroisonicotinic acid provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2083 mL	26.0417 mL	52.0833 mL
5 mM	1.0417 mL	5.2083 mL	10.4167 mL
10 mM	0.5208 mL	2.6042 mL	5.2083 mL
50 mM	0.1042 mL	0.5208 mL	1.0417 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481