


C.I. 20195

Chemical Properties

CAS No. : 5850-16-8
 Formula: C₂₆H₁₆N₄Na₂O₈S₂
 Molecular Weight: 624.55
 Storage: Keep away from direct sunlight
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	C.I. 20195 is a brown dye.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6012 mL	8.0058 mL	16.0115 mL
5 mM	0.3202 mL	1.6012 mL	3.2023 mL
10 mM	0.1601 mL	0.8006 mL	1.6012 mL
50 mM	0.032 mL	0.1601 mL	0.3202 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481