


Fipravorimat

Chemical Properties

CAS No. :	1818867-24-1
Formula:	C ₄₃ H ₆₇ N ₂ O ₄ S
Molecular Weight:	727.07
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fipravorimat is a potent inhibitor of HIV-1 with potential applications in HIV and AIDS research.
Targets(IC50)	Others,HIV Protease

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3754 mL	6.8769 mL	13.7538 mL
5 mM	0.2751 mL	1.3754 mL	2.7508 mL
10 mM	0.1375 mL	0.6877 mL	1.3754 mL
50 mM	0.0275 mL	0.1375 mL	0.2751 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481