

(S,S)-D211

Chemical Properties


CAS No. : 1971075-99-6

Formula: C₃₉H₃₆N₄O₆

Molecular Weight: 656.73

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S,S)-D211 is a DNA damaging payload.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5227 mL	7.6135 mL	15.227 mL
5 mM	0.3045 mL	1.5227 mL	3.0454 mL
10 mM	0.1523 mL	0.7613 mL	1.5227 mL
50 mM	0.0305 mL	0.1523 mL	0.3045 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481