

2-[[4-ethenyl-5-hydroxy-3-(2-hydroxyethyl)pentyl]oxy]-6-(hydroxymethyl)oxane-3,4,5-triol

Chemical Properties

CAS No. :

Formula: C15H28O8

Molecular Weight: 336.38

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(2R,3R,4S,5S,6R)-2-[[4-ethenyl-5-hydroxy-3-(2-hydroxyethyl)pentyl]oxy]-6-(hydroxymethyl)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9728 mL	14.8641 mL	29.7283 mL
5 mM	0.5946 mL	2.9728 mL	5.9457 mL
10 mM	0.2973 mL	1.4864 mL	2.9728 mL
50 mM	0.0595 mL	0.2973 mL	0.5946 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481