



## Piadar

## Chemical Properties

CAS No. :	78436-02-9
Formula:	C31H40O7
Molecular Weight:	524.654
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Piadar is a new antimycotic containing the antifungal agent siccanin.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.906 mL	9.5302 mL	19.0603 mL
5 mM	0.3812 mL	1.906 mL	3.8121 mL
10 mM	0.1906 mL	0.953 mL	1.906 mL
50 mM	0.0381 mL	0.1906 mL	0.3812 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Meyer-Rohn J. [Piadar (Ladar): new antimycotic drug containing the antifungal agent siccanin (author's transl)]. Mykosen. 1979 Jul;22(7):255-8. German. PubMed PMID: 381921.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481