


2',3',5'-Tri-O-benzoyl-5-methoxyuridine

Chemical Properties

CAS No. :	37805-86-0
Formula:	C ₃₁ H ₂₆ N ₂ O ₁₀
Molecular Weight:	586.55
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2',3',5'-Tri-O-benzoyl-5-methoxyuridine is a 5-Modified Pyrimidine Nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7049 mL	8.5244 mL	17.0488 mL
5 mM	0.341 mL	1.7049 mL	3.4098 mL
10 mM	0.1705 mL	0.8524 mL	1.7049 mL
50 mM	0.0341 mL	0.1705 mL	0.341 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481