

## (2-Phenyl-1,3-thiazol-4-yl)methanol

## Chemical Properties

CAS No. :	23780-13-4
Formula:	C10H9NOS
Molecular Weight:	191.25
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	(2-Phenyl-1,3-thiazol-4-yl)methanol ,with CAS No. 23780-13-4, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. (2-Phenyl-1,3-thiazol-4-yl)methanol provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2288 mL	26.1438 mL	52.2876 mL
5 mM	1.0458 mL	5.2288 mL	10.4575 mL
10 mM	0.5229 mL	2.6144 mL	5.2288 mL
50 mM	0.1046 mL	0.5229 mL	1.0458 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481