


## 1,2,3,4-Tetrahydro-6-quinolinecarboxylic Acid

## Chemical Properties

CAS No. :	5382-49-0
Formula:	C10H11NO2
Molecular Weight:	177.203
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1,2,3,4-Tetrahydro-6-quinolinecarboxylic Acid is a useful organic compound for research related to life sciences. The catalog number is Fr21307 and the CAS number is 5382-49-0.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: 35.44 mg/mL (200 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---------------------------------------------------------------------------------------------------------------------------

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.6433 mL	28.2167 mL	56.4334 mL
5 mM	1.1287 mL	5.6433 mL	11.2867 mL
10 mM	0.5643 mL	2.8217 mL	5.6433 mL
50 mM	0.1129 mL	0.5643 mL	1.1287 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481