

Compound N050-0020

Chemical Properties


CAS No. :

Formula: C₂₄H₃₆O₄

Molecular Weight: 388.548

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N050-0020 is a useful organic compound for research related to life sciences and the catalog number is T131248.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5737 mL	12.8684 mL	25.7367 mL
5 mM	0.5147 mL	2.5737 mL	5.1473 mL
10 mM	0.2574 mL	1.2868 mL	2.5737 mL
50 mM	0.0515 mL	0.2574 mL	0.5147 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481