

Erythrinin H

Chemical Properties

CAS No. :	1616592-62-1
Formula:	C17H12O7
Molecular Weight:	328.27
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Erythrinin H is a useful organic compound for research related to life sciences and the catalog number is T125680.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0463 mL	15.2314 mL	30.4627 mL
5 mM	0.6093 mL	3.0463 mL	6.0925 mL
10 mM	0.3046 mL	1.5231 mL	3.0463 mL
50 mM	0.0609 mL	0.3046 mL	0.6093 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481