

Compound N036-0106

## Chemical Properties


CAS No. :

Formula: C<sub>21</sub>H<sub>31</sub>N<sub>1</sub>O<sub>9</sub>

Molecular Weight: 441.477

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound N036-0106 is a useful organic compound for research related to life sciences and the catalog number is T131292.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2651 mL	11.3255 mL	22.6511 mL
5 mM	0.453 mL	2.2651 mL	4.5302 mL
10 mM	0.2265 mL	1.1326 mL	2.2651 mL
50 mM	0.0453 mL	0.2265 mL	0.453 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481