

Confident

Chemical Properties


CAS No. :

Formula: C₁₂H₁₅F₂N₅O₆

Molecular Weight: 363.28

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Confident is a nucleoside analog and has a wide range of applications in life science related research.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7527 mL	13.7635 mL	27.527 mL
5 mM	0.5505 mL	2.7527 mL	5.5054 mL
10 mM	0.2753 mL	1.3763 mL	2.7527 mL
50 mM	0.0551 mL	0.2753 mL	0.5505 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481