


Ritanserin tartrate

Chemical Properties

CAS No. :	93076-39-2
Formula:	C ₃₁ H ₃₁ F ₂ N ₃ O ₇ S
Molecular Weight:	627.66
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ritanserin tartrate is a selective and potent serotonin-2 antagonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5932 mL	7.9661 mL	15.9322 mL
5 mM	0.3186 mL	1.5932 mL	3.1864 mL
10 mM	0.1593 mL	0.7966 mL	1.5932 mL
50 mM	0.0319 mL	0.1593 mL	0.3186 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481