

tetradec-5-ynoic acid

Chemical Properties

CAS No. :

Formula: C₁₄H₂₄O₂

Molecular Weight: 224.34

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	tetradec-5-ynoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4575 mL	22.2876 mL	44.5752 mL
5 mM	0.8915 mL	4.4575 mL	8.915 mL
10 mM	0.4458 mL	2.2288 mL	4.4575 mL
50 mM	0.0892 mL	0.4458 mL	0.8915 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481