


Beiwutine

Chemical Properties

CAS No. :	76918-93-9
Formula:	C33H45NO12
Molecular Weight:	647.718
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Beiwutine is a natural product for research related to life sciences. The catalog number is TN6285 and the CAS number is 76918-93-9.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5439 mL	7.7194 mL	15.4388 mL
5 mM	0.3088 mL	1.5439 mL	3.0878 mL
10 mM	0.1544 mL	0.7719 mL	1.5439 mL
50 mM	0.0309 mL	0.1544 mL	0.3088 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481