


## Boc-NH-PEG5-azide

## Chemical Properties

CAS No. :	911209-07-9
Formula:	C17H34N4O7
Molecular Weight:	406.47
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Boc-NH-PEG5-azide, a PEG-based linker for PROTACs, connects two essential ligands, facilitating the formation of PROTAC molecules and enabling selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs utilize two distinct ligands connected by a linker, one binding to an E3 ubiquitin ligase and the other to the target protein, thereby harnessing the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4602 mL	12.301 mL	24.6021 mL
5 mM	0.492 mL	2.4602 mL	4.9204 mL
10 mM	0.246 mL	1.2301 mL	2.4602 mL
50 mM	0.0492 mL	0.246 mL	0.492 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481