


JMV3008

## Chemical Properties

CAS No. : 925239-09-4  
 Formula: C<sub>36</sub>H<sub>39</sub>N<sub>7</sub>O<sub>3</sub>  
 Molecular Weight: 617.74  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	JMV3008 is a bioactive chemical.
Targets(IC50)	Others,GHSR

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6188 mL	8.094 mL	16.188 mL
5 mM	0.3238 mL	1.6188 mL	3.2376 mL
10 mM	0.1619 mL	0.8094 mL	1.6188 mL
50 mM	0.0324 mL	0.1619 mL	0.3238 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481