

## Litsomentol

## Chemical Properties

CAS No. :

Formula: C<sub>32</sub>H<sub>54</sub>O<sub>3</sub>

Molecular Weight: 486.781

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Litsomentol is a useful organic compound for research related to life sciences and the catalog number is T131514.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0543 mL	10.2716 mL	20.5432 mL
5 mM	0.4109 mL	2.0543 mL	4.1086 mL
10 mM	0.2054 mL	1.0272 mL	2.0543 mL
50 mM	0.0411 mL	0.2054 mL	0.4109 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481