


## Tetrahydropalmatine

## Chemical Properties

CAS No. :	2934-97-6
Formula:	C <sub>21</sub> H <sub>25</sub> NO <sub>4</sub>
Molecular Weight:	355.43
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Tetrahydropalmatine (DL-Tetrahydropalmatine) is an active component isolated from corydalis, suppressing amygdaloid release of dopamine to inhibit an epileptic attack in rats.
Targets(IC50)	Apoptosis,Dopamine Receptor

## Solubility Information

Solubility	DMSO: 16.7 mg/mL (46.99 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+90% Saline: 1.67 mg/mL (4.7 mM),Solution. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8135 mL	14.0675 mL	28.1349 mL
5 mM	0.5627 mL	2.8135 mL	5.627 mL
10 mM	0.2813 mL	1.4067 mL	2.8135 mL
50 mM	0.0563 mL	0.2813 mL	0.5627 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

XING Shu-Hua, et al. Acta Pharmacologica Sinica, 1994-01

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481