

Compound T128327

Chemical Properties


CAS No. :

Formula: C₂₉H₄₀N₂O₅

Molecular Weight: 496.65

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-[(6,7-dimethoxy-1,2,3,4-tetrahydroisoquinolin-1-yl)methyl]-3-ethyl-9,10-dimethoxy-1H,2H,3H,4H,5H,6H,7H,11bH-pyrido[2,1-a]isoquinolin-5-ium-5-olate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0135 mL	10.0675 mL	20.1349 mL
5 mM	0.4027 mL	2.0135 mL	4.027 mL
10 mM	0.2013 mL	1.0067 mL	2.0135 mL
50 mM	0.0403 mL	0.2013 mL	0.4027 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481