


Conduritol A

Chemical Properties

CAS No. :	526-87-4
Formula:	C6H10O4
Molecular Weight:	146.14
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Conduritol A has a hypoglycemic effect, can have an effect on regulating the metabolism of blood lipid, free-radical scavenging, enhancing the antioxidant ability, potentiating immune function; it also can markedly prevent the diabetic rats from getting cataracts through the inhibition of lens aldose reductase.
Targets(IC50)	TNF

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.8428 mL	34.2138 mL	68.4275 mL
5 mM	1.3686 mL	6.8428 mL	13.6855 mL
10 mM	0.6843 mL	3.4214 mL	6.8428 mL
50 mM	0.1369 mL	0.6843 mL	1.3686 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481