

Peniprequinolone

Chemical Properties

CAS No. :

Formula: C₂₂H₂₅N₅O₅

Molecular Weight: 383.444

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Peniprequinolone is a useful organic compound for research related to life sciences and the catalog number is T130829.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.608 mL	13.0398 mL	26.0797 mL
5 mM	0.5216 mL	2.608 mL	5.2159 mL
10 mM	0.2608 mL	1.304 mL	2.608 mL
50 mM	0.0522 mL	0.2608 mL	0.5216 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481