

exo BCN-O-PNB

Chemical Properties

CAS No. : 1380006-72-3


Formula: C17H17NO5

Molecular Weight: 315.32

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	exo BCN-O-PNB is a alkyl/ether-based linker for PROTACs which joins two essential ligands, crucial for forming PROTAC molecules. This linker enables selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs utilize the intracellular ubiquitin-proteasome system for the selective degradation of target proteins by linking two distinct ligands: one for an E3 ubiquitin ligase and one for the target protein[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1714 mL	15.8569 mL	31.7138 mL
5 mM	0.6343 mL	3.1714 mL	6.3428 mL
10 mM	0.3171 mL	1.5857 mL	3.1714 mL
50 mM	0.0634 mL	0.3171 mL	0.6343 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481