


ST-115

## Chemical Properties

CAS No. : 1840873-90-6  
 Formula: C<sub>21</sub>H<sub>32</sub>FN<sub>3</sub>O<sub>5</sub>  
 Molecular Weight: 457.56  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	ST-115 is a potent and specific inhibitor of aminopeptidase P2. When ST-115 is administered intravenously at the start of reperfusion, it reduces bradykinin degradation. This increases bradykinin's concentration in the capillaries and enhances its protective effects.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1855 mL	10.9275 mL	21.8551 mL
5 mM	0.4371 mL	2.1855 mL	4.371 mL
10 mM	0.2186 mL	1.0928 mL	2.1855 mL
50 mM	0.0437 mL	0.2186 mL	0.4371 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481