


Suberosol

Chemical Properties

CAS No. :	151368-42-2
Formula:	C ₃₁ H ₅₀ O ₂
Molecular Weight:	454.73
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Suberosol possesses anti-HIV replication activity.
Targets(IC50)	Others,HIV Protease

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1991 mL	10.9955 mL	21.9911 mL
5 mM	0.4398 mL	2.1991 mL	4.3982 mL
10 mM	0.2199 mL	1.0996 mL	2.1991 mL
50 mM	0.044 mL	0.2199 mL	0.4398 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Li HY, Sun NJ, Kashiwada Y, Sun L, Snider JV, Cosentino LM, Lee KH. Anti-AIDS agents, 9. Suberosol, a new C₃₁ lanostane-type triterpene and anti-HIV principle from *Polyalthia suberosa*. *J Nat Prod*. 1993 Jul;56(7):1130-3. PubMed PMID: 8377018.

Zhao W, Fan C, Qin G. [A study on major chemical components of *Uvaria grandiflora*]. *Zhongguo Zhong Yao Za Zhi*. 1999 Aug;24(8):476-7, 510-1. Chinese. PubMed PMID: 12205865.

Yang PS, Cheng MJ, Peng CF, Chen JJ, Chen IS. Endiandric Acid Analogues from the Roots of *Beilschmiedia erythrophloia*. *J Nat Prod*. 2009 Jan;72(1):53-8. doi: 10.1021/np800504w. PubMed PMID: 19072217.

Zhang W, Tang SA, Duan HQ. [Studies on chemical constituents of rhizome of *Buddleia davidii*]. *Zhong Yao Cai*. 2009 Apr;32(4):515-7. Chinese. PubMed PMID: 19645234.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481