


## Axillaridine A

## Chemical Properties

CAS No. :	128255-16-3
Formula:	C30H42N2O2
Molecular Weight:	462.67
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Axillaridine A is a new cholinesterase inhibitors, it may act as potential leads in the discovery of clinically useful inhibitors for nervous-system disorders, particularly by reducing memory deficiency in Alzheimer's disease patients by potentiating and effecting the cholinergic transmission process. (+)-Axillaridine A has significant activity as antiestrogen binding site (AEBS)-inhibitory agents.
Targets(IC50)	Others,Cholinesterase (ChE)

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1614 mL	10.8068 mL	21.6137 mL
5 mM	0.4323 mL	2.1614 mL	4.3227 mL
10 mM	0.2161 mL	1.0807 mL	2.1614 mL
50 mM	0.0432 mL	0.2161 mL	0.4323 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481