

Paludosine

Chemical Properties

CAS No. :

Formula: C₃₁H₄₅N₀₈

Molecular Weight: 559.7

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Paludosine is a natural product that can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7867 mL	8.9334 mL	17.8667 mL
5 mM	0.3573 mL	1.7867 mL	3.5733 mL
10 mM	0.1787 mL	0.8933 mL	1.7867 mL
50 mM	0.0357 mL	0.1787 mL	0.3573 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481