


4-O-D-Glucopyranosyl-p-coumaric acid

Chemical Properties

CAS No. :	117405-49-9
Formula:	C ₁₅ H ₁₈ O ₈
Molecular Weight:	326.301
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-O-D-Glucopyranosyl-p-coumaric acid is a useful organic compound for research related to life sciences. The catalog number is T126288 and the CAS number is 117405-49-9.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0647 mL	15.3233 mL	30.6466 mL
5 mM	0.6129 mL	3.0647 mL	6.1293 mL
10 mM	0.3065 mL	1.5323 mL	3.0647 mL
50 mM	0.0613 mL	0.3065 mL	0.6129 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481