


ISAM-140

## Chemical Properties

CAS No. : 932191-62-3  
 Formula: C<sub>19</sub>H<sub>19</sub>N<sub>3</sub>O<sub>3</sub>  
 Molecular Weight: 337.37  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	ISAM-140 is a potent and highly selective antagonist of A2B adenosine receptor.
Targets(IC50)	Adenosine Receptor

## Solubility Information

Solubility	DMSO: 3.38 mg/mL (10.02 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9641 mL	14.8205 mL	29.641 mL
5 mM	0.5928 mL	2.9641 mL	5.9282 mL
10 mM	0.2964 mL	1.4821 mL	2.9641 mL
50 mM	0.0593 mL	0.2964 mL	0.5928 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

El Maatougui A, et al. Discovery of Potent and Highly Selective A2B Adenosine Receptor Antagonist Chemotypes. J Med Chem. 2016 Mar 10;59(5):1967-83.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481