

## 4-(2-Nitroethyl)phenyl primeveroside

## Chemical Properties


CAS No. :

Formula: C<sub>19</sub>H<sub>27</sub>N<sub>2</sub>O<sub>12</sub>

Molecular Weight: 461.42

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-(2-Nitroethyl)phenyl primeveroside is a useful organic compound for research related to life sciences and the catalog number is T130388.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1672 mL	10.8361 mL	21.6722 mL
5 mM	0.4334 mL	2.1672 mL	4.3344 mL
10 mM	0.2167 mL	1.0836 mL	2.1672 mL
50 mM	0.0433 mL	0.2167 mL	0.4334 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481