


Acaprazine

Chemical Properties

CAS No. : 55485-20-6
 Formula: C₁₅H₂₁Cl₂N₃O
 Molecular Weight: 330.25
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Acaprazine is an adrenergic and anxiolytic drug in the class of phenylpiperazine.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.028 mL	15.140 mL	30.2801 mL
5 mM	0.6056 mL	3.028 mL	6.056 mL
10 mM	0.3028 mL	1.514 mL	3.028 mL
50 mM	0.0606 mL	0.3028 mL	0.6056 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481