


3,6-Dimethoxyapigenin

Chemical Properties

CAS No. :	22697-65-0
Formula:	C ₁₇ H ₁₄ O ₇
Molecular Weight:	330.29
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3,6-Dimethoxyapigenin is a natural product of Dodonaea, Sapindaceae. The catalog number is TN2907 and the CAS number is 22697-65-0. 3,6-Dimethoxyapigenin can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0276 mL	15.1382 mL	30.2764 mL
5 mM	0.6055 mL	3.0276 mL	6.0553 mL
10 mM	0.3028 mL	1.5138 mL	3.0276 mL
50 mM	0.0606 mL	0.3028 mL	0.6055 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481