


TMC 120C

Chemical Properties

CAS No. : 250231-82-4
 Formula: C₁₅H₁₅NO₃
 Molecular Weight: 257.289
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	TMC 120C is a natural product that can be used as a reference standard. The CAS number of TMC 120C is 250231-82-4.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8867 mL	19.4333 mL	38.8666 mL
5 mM	0.7773 mL	3.8867 mL	7.7733 mL
10 mM	0.3887 mL	1.9433 mL	3.8867 mL
50 mM	0.0777 mL	0.3887 mL	0.7773 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481