

## Compound T128487

## Chemical Properties


CAS No. :

Formula: C35H56O8

Molecular Weight: 604.83

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(2S,4aR,6bR,10S,12aR,14bS)-2-(hydroxymethyl)-2,6a,6b,9,9,12a-hexamethyl-10-(((2S,3R,4S,5S)-3,4,5-trihydroxyoxan-2-yl]oxy)-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,12b,13,14b-icosahydricene-4a-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6534 mL	8.2668 mL	16.5336 mL
5 mM	0.3307 mL	1.6534 mL	3.3067 mL
10 mM	0.1653 mL	0.8267 mL	1.6534 mL
50 mM	0.0331 mL	0.1653 mL	0.3307 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481