


Sonopiprazole hydrochloride

Chemical Properties

CAS No. :	170857-36-0
Formula:	C ₂₁ H ₂₈ ClN ₃ O ₃ S
Molecular Weight:	437.98
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sonopiprazole HCl is a selective dopamine D4 antagonist.
Targets(IC50)	Others, Dopamine Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2832 mL	11.416 mL	22.8321 mL
5 mM	0.4566 mL	2.2832 mL	4.5664 mL
10 mM	0.2283 mL	1.1416 mL	2.2832 mL
50 mM	0.0457 mL	0.2283 mL	0.4566 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481