

SPDP-PEG12-NHS ester

Chemical Properties

CAS No. :

Formula: C39H65N3O17S2

Molecular Weight: 912.07



Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SPDP-PEG12-NHS ester is a PEG-based linker for PROTACs that joins two essential ligands, crucial for forming PROTAC molecules, enabling selective protein degradation via the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands linked together; one binds to an E3 ubiquitin ligase and the other to a target protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0964 mL	5.482 mL	10.9641 mL
5 mM	0.2193 mL	1.0964 mL	2.1928 mL
10 mM	0.1096 mL	0.5482 mL	1.0964 mL
50 mM	0.0219 mL	0.1096 mL	0.2193 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481