

CBR Agonist-1

Chemical Properties

CAS No. :

Formula: C27H27FN4O

Molecular Weight: 442.53

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	CBR Agonist-1 (27a-cis) is a cannabinoid receptor (CBR) agonist that targets CB1R with a K_i of 0.18 μM and CB2R with a K_i of 1.22 μM . It can be utilized to investigate diseases related to the endogenous cannabinoid system.
Targets(IC50)	Cannabinoid Receptor,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2597 mL	11.2987 mL	22.5973 mL
5 mM	0.4519 mL	2.2597 mL	4.5195 mL
10 mM	0.226 mL	1.1299 mL	2.2597 mL
50 mM	0.0452 mL	0.226 mL	0.4519 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481