


Propargyl p-toluenesulfonate

Chemical Properties

CAS No. :	6165-76-0
Formula:	C ₁₀ H ₁₀ O ₃ S
Molecular Weight:	210.25
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Propargyl p-toluenesulfonate is a useful organic compound for research related to life sciences. The catalog number is T124543 and the CAS number is 6165-76-0.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.7562 mL	23.7812 mL	47.5624 mL
5 mM	0.9512 mL	4.7562 mL	9.5125 mL
10 mM	0.4756 mL	2.3781 mL	4.7562 mL
50 mM	0.0951 mL	0.4756 mL	0.9512 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481