


PPAR γ agonist 7

Chemical Properties

CAS No. :	2569295-93-6
Formula:	C ₂₀ H ₃₀ O ₆
Molecular Weight:	366.45
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	PPAR γ agonist 7 (Compound 3a) is a highly potent and selective peroxisome proliferator-activated receptor gamma (PPAR γ) agonist, which specifically stimulates adiponectin production in human bone marrow mesenchymal stem cells (hBM-MSCs), demonstrating an innovative full agonist profile with an EC ₅₀ value of 4.34 μ M [1].
Targets(IC50)	Others,PPAR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7289 mL	13.6444 mL	27.2889 mL
5 mM	0.5458 mL	2.7289 mL	5.4578 mL
10 mM	0.2729 mL	1.3644 mL	2.7289 mL
50 mM	0.0546 mL	0.2729 mL	0.5458 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481