


## Butyl carbitol acetate

## Chemical Properties

CAS No. :	124-17-4
Formula:	C10H20O4
Molecular Weight:	204.27
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Butyl carbitol acetate enables deep permeation into the textile, with its low vapor pressure and boiling point.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8955 mL	24.4774 mL	48.9548 mL
5 mM	0.9791 mL	4.8955 mL	9.791 mL
10 mM	0.4895 mL	2.4477 mL	4.8955 mL
50 mM	0.0979 mL	0.4895 mL	0.9791 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481