


C.I. Solvent Green 3

Chemical Properties

CAS No. :	128-80-3
Formula:	C ₂₈ H ₂₂ N ₂ O ₂
Molecular Weight:	418.49
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	C.I. Solvent Green 3 is a green dye, an is a green dye, an anthraquinone derivative..
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3895 mL	11.9477 mL	23.8954 mL
5 mM	0.4779 mL	2.3895 mL	4.7791 mL
10 mM	0.239 mL	1.1948 mL	2.3895 mL
50 mM	0.0478 mL	0.239 mL	0.4779 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481