


## Compound FL0045

## Chemical Properties

CAS No. :	1628184-67-7
Formula:	C <sub>5</sub> H <sub>6</sub> F <sub>3</sub> NO
Molecular Weight:	153.1
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Compound FL0045, with CAS No. 1628184-67-7, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Compound FL0045 provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.5317 mL	32.6584 mL	65.3168 mL
5 mM	1.3063 mL	6.5317 mL	13.0634 mL
10 mM	0.6532 mL	3.2658 mL	6.5317 mL
50 mM	0.1306 mL	0.6532 mL	1.3063 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481