

(11E,15Z)-9,10,13-trihydroxyoctadeca-11,15-dienoic acid

Chemical Properties

CAS No. :

Formula: C₁₈H₃₂O₅

Molecular Weight: 328.45

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(11E,15Z)-9,10,13-trihydroxyoctadeca-11,15-dienoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0446 mL	15.223 mL	30.446 mL
5 mM	0.6089 mL	3.0446 mL	6.0892 mL
10 mM	0.3045 mL	1.5223 mL	3.0446 mL
50 mM	0.0609 mL	0.3045 mL	0.6089 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481