

PPC-NB

Chemical Properties

CAS No. : 1887040-81-4


Formula: C₁₅H₁₄N₂O₅S₂

Molecular Weight: 366.41

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	PPC-NB is a glutathione-cleavable linker used in antibody-drug conjugates (ADC) [1].
Targets(IC50)	ADC Linker,PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7292 mL	13.6459 mL	27.2918 mL
5 mM	0.5458 mL	2.7292 mL	5.4584 mL
10 mM	0.2729 mL	1.3646 mL	2.7292 mL
50 mM	0.0546 mL	0.2729 mL	0.5458 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Robert A. Volkmann, et al. Compounds, compositions, and methods for treatment of diseases involving acidic or hypoxic diseased tissues. WO2019136298A1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481