


## 4-Methylsulfonylaniline

## Chemical Properties

CAS No. :	5470-49-5
Formula:	C7H9NO2S
Molecular Weight:	171.21
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-Methylsulfonylaniline ,with CAS No. 5470-49-5, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 4-Methylsulfonylaniline provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8408 mL	29.2039 mL	58.4078 mL
5 mM	1.1682 mL	5.8408 mL	11.6816 mL
10 mM	0.5841 mL	2.9204 mL	5.8408 mL
50 mM	0.1168 mL	0.5841 mL	1.1682 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481